


Parameter	Specification		
Description	Thin polyester heat sealing film, easily pierceable with auto sampler needles/ABI® 3730; suitable for PCR, qPCR, and optical applications.		
Composition	Thin clear polyester backed film that seals to PP, PS, PE & COC plates. Not suitable for peeling.		
Seal Type	Heat.		
Ordering Information	4ti-0580	610 m x 78 mm roll Compatible with Automated Roll Heat Sealer, Thermo Fisher ALPS 300™ / ALPS 3000™ / KBiosystemsWasp™, KBioscience FlexiSeal and Cube	
	4ti-0580/122	122 m x 78 mm roll Compatible with Automated Roll Heat Sealer with optional dust cover	
	4ti-0582	500 m x 115 mm roll Compatible with Agilent (Velocity 11) PlateLoc®	
	4ti-0581	Roll of 1000 perforated sheets (125 x 78 mm)	
Application	QPCR and other optical applications such as fluorescence or colorimetric measurements. The thin nature of this seal enables it to be easily pierced with an auto sampler needle or the probe from the ABI 3730 sequencer.		
Certification	Free from human genomic DNA, nucleases, and pyrogens.		
Seal Integrity Range	-20 °C to 80 °C or 110 °C with pressurised PCR heated lid.		
Storage Conditions	Room temperature (10-27 °C). Avoid direct exposure to light and high humidity.		
Material Safety Data	Non-hazardous		
Schematic Overview	<div style="text-align: center;">  <p>Outer Surface</p> <p>Sealing Surface</p> <p>Polyester 25 µm</p> <p>Copolymer 50 µm</p> </div>		
Pierce Force	~6.0 N (ASTM D1000)		
Sealing Parameters [^]	96-well PCR Plate	165-175 °C, 3 seconds	
	384-well PCR Plate	165-175 °C, 2 seconds	
	Vision Plate™	175-185 °C, 2-3 seconds	
Physical Properties	Thickness	Polyester Copolymer	25.00 µm 50.00 µm
	Seal Strength	>60 N/25 mm	

**These are recommended sealing conditions for the 4s3™ heat sealer.*

Proprietary information - This document and the information disclosed herein is confidential and proprietary to Azenta, Inc. and may not be reproduced in whole or in part or disclosed to any third party or used without the prior written consent of Azenta, Inc.