

Title: Impression Update	
Document No: NOC-2003	
Revision: 1	March 2020

We are writing to inform you that as part of our continuous improvement, the Impression MKI units (part number 20-2101-A) will now include a Linear Barcode Module as standard. This will be effective from Serial Number 2101-A-0120-001. The part number remains the same. The Impression MKI model 20-2101 and 20-2101-B will be made obsolete as part of this notice of change.

The linear barcode reader (70-2010) will still be available as a separate part.

1. Affected Products

Part number	Description	Update
20-2101-A	Impression Rapid Rack Scanner (MKI) with Opticon Linear Barcode Reader Small Form Factor Whole rack reader for racks of 48 / 96 2D labelled tubes; small-footprint scanner-based solution for rapid reading and integration. Small form factor scanner has a lower depth of focus and is not suitable for all racks	Now includes a Linear Barcode Reader (70-2010) as standard
20-2101	Impression Rapid Rack Scanner (MKI) Small Form Factor Whole rack reader for racks of 48 / 96 2D labelled tubes; small-footprint scanner-based solution for rapid reading and integration. Small form factor scanner has a lower depth of focus and is not suitable for all racks	Obsolete
20-2101-В	Impression Rapid Rack Scanner (MKI) with Opticon Linear Barcode Reader Small Form Factor Whole rack reader for racks of 48 / 96 2D labelled tubes; small-footprint scanner-based solution for rapid reading and integration. Small form factor scanner has a lower depth of focus and is not suitable for all racks	Obsolete
70-2010	(Optional) 1D Linear Barcode Module for Impression USB Linear Barcode Module for reading 1D barcodes on the side of an SBS rack. 1 per case	Still available





Figure 1: Impression MKI with Opticon Barcode Reader (part number 20-2101-A)

2. Reason for Change

Due to the MKI small form factor, it is not always possible to read 2D rack IDs on the base of SBS racks, therefore we have decided to include a linear barcode reader as part of the Impression MKI rapid rack scanner.

Part 20-2101-B and 20-2101 will be made obsolete as they are duplicate part numbers and unnecessary to maintain as part of our portfolio of products.

3. Date of Change

Effective 5th March 2020

Units with the serial number 2101-A-0120-001 onwards will include the Linear Barcode Reader.

If you require any further information please feel free to contact your local rep, distributor or email the office on <u>blss.sales@brooks.com</u>.