

Description	96 Well Assay Plate, Pyramid Flat Bottom
Creation Date	02.09.2020
Product Change Notice Number	NOC 2015
Revision	01
This Product Change Notice supersedes the following documents:	

1. Description and Reason for Change

We are writing to inform you that 96 Well Assay Plate, Pyramid Flat Bottom will be discontinued on the 2nd September 2020. The discontinuation of the product is in line with our practice of rationalizing low demand products and consolidating products with similar features.

96 Well Assay Plate can be used as a replacement. For further information about replacement product, please follow the below link -

<https://www.brookslifesciences.com/products/96-well-assay-plate>

2. Communication Distribution – Need to know list

Customer's Procurement	C&I Sales to inform
Customer's End User	C&I Sales to inform
Customer's Service Team	C&I Sales to inform
Brooks Life Sciences Service Team	Product Management to inform
Brooks Life Sciences Customer Care Team	Product Management to inform

3. Classification of Change

Customer Notification	Yes
Information Only	

4. Impact of Change

If you switch to alternative 96 Well Assay Plate, please note the following:

- There are some differences with shape of wells (replacement plates have round wells and standard flat bottom instead of square wells and pyramid flat bottom).
- Replacement plates comes in black (4ti-0263) and white (4ti-0273) colour options. We don't have replacement product for clear plate version.
- Replacement plates are made from polystyrene which is same polymer material that is used to produce 96 Well Assay Plate, Pyramid Flat Bottom.

5. Affected Part Numbers

Item Number	Description
4ti-0255	96 Well Assay Plate, Pyramid Flat Bottom; 1.1 ml square wells, clear PS, cut corner H12, 21 plates
4ti-0265	96 Well Assay Plate, Pyramid Flat Bottom; 1.1 ml square wells, black PS, cut corner H12, 21 plates
4ti-0275	96 Well Assay Plate, Pyramid Flat Bottom; 1.1 ml square wells, white PS, cut corner H12, 21 plates

Please note this affects all custom and treated products which contain any of these base part numbers.

6. Replacement Products

Item Number	Alternative Item Number	Alternative Description
4ti-0255	N/A	No suitable replacement available
4ti-0265	4ti-0263	96 Well Assay Plate; 0.35 ml round wells, flat base, black PS, cut corner A1/H1, 100 plates
4ti-0275	4ti-0273	96 Well Assay Plate; 0.35 ml round wells, flat base, white PS, cut corner A1/H1, 100 plates

7. Recommended Action

Limited stock of 96 Well Assay Plate, Pyramid Flat Bottom is available at 50% discount while stocks last.

Further information about replacement product, 96 Well Assay Plate is available on our website. Please follow the link below for technical drawings, specifications and sample request –

<https://www.brookslifesciences.com/products/96-well-assay-plate>

8. Specification of replacement products

Parameter	96 Well Assay Plate, Pyramid Flat Bottom	96 Well Assay Plate
Plate length	127.76 ± 0.25mm	127.76 ± 0.25mm
Plate width	85.45 ± 0.25mm	85.48 ± 0.25mm
Plate height	19.60 ± 0.25mm	14.35 ± 0.25mm
Well depth	17.50 ± 0.25mm	10.8 ± 0.25mm
Well diameter/length	8.2 ± 0.10mm	6.3 ± 0.10mm
Distance to center of A1 from top edge	11.22 ± 0.25mm	11.24 ± 0.25mm
Distance to center of A1 from left edge	14.38 ± 0.25mm	14.38 ± 0.25mm
Pitch (distance between A1 and A2)	9.00mm	9.00mm
Well Bottom	Pyramid flat bottom	Flat bottom
Top Well shape	Square	Round

Table 8.1. comparison of affected and replacement products

9. Timeframe for implementation

Effective from 2nd September 2020

10. For more information please contact Brooks Life Sciences

If you require further information, please contact your local Account or Distribution Manager, or call our Customer Service team on

North America +1.800.379.7221
Europe +44.0.161.777.2000
Japan +81.45.4477.5570 (ext. 24)