


DOCUMENT NUMBER: QMS101778	TITLE: Notice of Change Form	
REVISION: A DCO# DCO005185	DOCUMENT CLASSIFICATION: 04-Form, Template or Other	

Product Change Notice

Description	4ti-0750/TA/12 tooling updates
Creation Date	24th June 2022
Product Change Notice Number	2208
Revision	Revision A
This Product Change Notice supersedes the following documents:	N/A

1. Description and Reason for Change

We are writing to inform you that as part of our aim for continuous improvement and to best facilitate customer demand, there will be minor changes to our 4ti-0750/TA/12 plates due to new tooling.

The primary difference is the deck of the plate will be smaller and as such, the alphanumerical grid references will be removed. As this product is used in strip/segment format – we have intentionally removed these references to reduce the amount of unnecessary plastic in the plate, ultimately resulting in less wastage. Offering a more environmentally friendly product that simultaneously eliminates manufacturing difficulties posed by previous tooling.

Additionally, the new design requires for a user to ‘tear’ 12 well strips away from the main body of the plate, rather than ‘break’ or ‘snap’ 12 wells strips off. The original design occasionally led to issues such as uneven breakage, or potential damage to other parts of the frame, deck, or wells of the plates. This new design mitigates against these issues, alongside offering a greater level of useability and reproducibility.


This change has become necessary as it is becoming increasingly difficult to maintain the quality and productivity levels of this SKU due to ageing production line equipment.

The newly introduced tool will be effective immediately, with details of the changes outlined in **Section 8**. Samples are available for evaluation to existing customers using the 4ti-0750/TA/12 plates.

2. Communication Distribution – Need to know list

Customer’s Procurement	Sales Managers to inform
Customer’s End User	Sales Managers to inform
Customer’s Service Team	Sales Managers to inform
Azenta Life Sciences Service Team	Product Management to inform
Azenta Life Sciences Customer Care Team	Product Management to inform

Confidential: The information is confidential and is to be used only in connection with matters authorized by Azenta and no part of it is to be disclosed to others without prior written permission from Azenta.

DOCUMENT NUMBER: QMS101778	TITLE: Notice of Change Form	
REVISION: A DCO# DCO005185	DOCUMENT CLASSIFICATION: 04-Form, Template or Other	

3. Classification of Change

Customer Notification	Yes
Information Only	

4. Impact of Change

As a result of tooling changes, there will be minor changes to the plates appearance. The deck of the plate will be smaller, and the alphanumeric references will be removed with no expected impact to intended use or compatibility.

5. Affected Part Numbers

Item #	Description
4ti-0750/TA/12	96 Well Non-Skirted PCR Plate, clear PP, breakable horizontally, high profile, cut corner H12, 50 plates per case.

Please note this affects all custom and treated products that contain this as a base part number.

6. Replacement Products

Not applicable.

7. Recommended Action


To review the differences detailed in section 8. We are happy to support and facilitate sample products of the preceding design for evaluation purposes if desired. These can be requested through the contact information listed in section 10.

Confidential: The information is confidential and is to be used only in connection with matters authorized by Azenta and no part of it is to be disclosed to others without prior written permission from Azenta.

DOCUMENT NUMBER: QMS101778	TITLE: Notice of Change Form	
REVISION: A DCO# DCO005185	DOCUMENT CLASSIFICATION: 04-Form, Template or Other	

8. Specification of current & replacement products


PREVIOUS DESIGN:



Confidential: The information is confidential and is to be used only in connection with matters authorized by Azenta and no part of it is to be disclosed to others without prior written permission from Azenta.

DOCUMENT NUMBER: QMS101778	TITLE: Notice of Change Form	
REVISION: A DCO# DCO005185	DOCUMENT CLASSIFICATION: O4-Form, Template or Other	


NEW DESIGN:



9. Timeframe for implementation

The new design will be implemented immediately.

Confidential: The information is confidential and is to be used only in connection with matters authorized by Azenta and no part of it is to be disclosed to others without prior written permission from Azenta.

DOCUMENT NUMBER: QMS101778	TITLE: Notice of Change Form	
REVISION: A DCO# DC0005185	DOCUMENT CLASSIFICATION: 04-Form, Template or Other	

10. For more information, please contact Azenta Life Sciences

If you require further information, please contact your local Account or Distribution Manager, or call our Customer Service team on:

North America +1.800.379.7221
Europe +44.0.161.777.2000
Japan +81.45.4477.5570 (ext. 24)

Confidential: The information is confidential and is to be used only in connection with matters authorized by Azenta and no part of it is to be disclosed to others without prior written permission from Azenta.