

CryoPod LN2 Filling Station Support Plans

Azenta Life Sciences products are designed and tested to suit the workflow demands of our customers. We know that some customers operate our instruments in environments where even the smallest downtime can be costly.

In recognition of this, Azenta Life Sciences offers standard Support Plans for the CryoPod® LN2 Filling Station to meet your workflow needs. These plans have been structured to offer flexibility so you can choose a package most appropriate to you. There are a range of benefits to offer peace of mind.



Key Benefits

- Mitigate downtime in your workflows
- Direct access to expert technical support
- Confidence in the protection of workflow
- Discounts for multi-year, multi-instrument contracts
- We can also provide tailored service options for customers with a more bespoke requirement

Validating the Installation

Your filling station should undergo a set of qualification procedures to ensure that it is installed correctly and meets all specifications:

- Installation Qualification (IQ): Oversight and verification of every physical aspect of the equipment (materials, dimensions, pressure ratings) and components (operational parameters, accuracy, voltage)
- Operational Qualification (OQ): Testing each individual component, feature, and physical specification of the equipment
- Performance Qualification (PQ) or IQ/OQ/PQ: Testing the operational requirements of the equipment under real-world conditions



AZENTA
LIFE SCIENCES

azenta.com



CryoPod LN2 Filling Station Support Plans

Support Plans

	Factory Warranty	Basic Package	Performance Package
CORRECTIVE SERVICES			
Unlimited phone & email support	✓	✓	✓
Priority on-site response commitment	Standard response time	Standard response time	3 business days
On-site corrective maintenance	✓	10% discount	✓
Warranty	1 Year: includes onsite service, travel, and spare parts		
REMOTE DIAGNOSTIC & REPAIR			
Unlimited priority phone & email support	✓	✓	✓
Remote Access Support, where possible	✓	✓	✓
PREVENTATIVE MAINTENANCE SERVICES (GASKET REPLACEMENT INCLUDED)			
Pre-scheduled, on-site preventive maintenance (PM) & PM parts kit (1/system/year)		✓	✓
Software & firmware updates during preventive maintenance, as required			✓
ADDITIONAL SERVICES AVAILABLE FOR PURCHASE			
Loaner equipment, based on availability	Not available	10% discount	15% discount
Loyalty discount on parts & accessories	Not available	10% discount	Included
Discount of on-site or institute-based training, where applicable	Full charges apply	5% discount	10% discount
Qualification & calibration services	Full charges apply	5% discount	10% discount

Part Numbers

- IOPQ documents | Part # 340594
- Service to perform the IOPQ | Part # 340592
- Gasket replacement included with preventative maintenance | Part # 308641



AZENTA
LIFE SCIENCES

azenta.com

© 2025 Azenta US, Inc. All rights reserved. All trademarks are property of Azenta US, Inc. unless otherwise specified. 25002-SD-2 0725