

NOT TO SCALE

127.76 ±0.25

103.50 ±0.25

12.13  
±0.25

4.50 typ

14.35  
±0.25

67.50 ±0.25


4.50 typ

8.99  
±0.25

X

X


85.48 ±0.25



2.50


3.70 typ

11.35



SECTION X-X

This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.



**AZENTA**  
LIFE SCIENCES

www.azenta.com

TITLE:

VISION PLATE 384

DATE:

22/3/07

SHEET No.:

1 OF 1

ISSUE:

A

All dimensions in millimetres