



NOTE:

1. SURFACE FINISH TO BE SMOOTH HIGH GLOSS - MIRROR POLISHED ON TOOL BY VDI 1
2. THE WHOLE PART MUST BE FLAT AND PARALLEL TO WITHIN ± 0.15
3. MAXIMUM PERMISSIBLE SIDE WALL BOWING ± 0.15
4. THE UNDIMENSIONED VALUES BY 3D CAD MODEL

ISSUE	DATE	AMENDMENT



 DO NOT SCALE IF IN DOUBT ASK	CUSTOMER	DRAWN	ZE	TITLE Standard PCR plate, semi-skirted with upstand MARKETING DRAWING
	CUSTOMER DRG No.	ISSUE	DATE	
	MATERIAL	POLYPROPYLENE Purell HP371P		CHECKED & SIGNATURE
	TOLERANCES	0 ±0.50 0.0 ±0.10 0.00 ±0.05 ANGULAR ±0°30'		DATE
	GENERAL TOLERANCES (UNLESS OTHERWISE STATED)	(ALL DIMENSIONS IN MILLIMETRES)		SCALE
		2:1		DRAWING No. 4ti-0735
				SHEET No. 1 OF 1
				ISSUE B